

## Assay Performance Data

<b>Naam assay:</b>	Thyroïdperoxidase IgG As [eenheden/volume] in serum of plasma [LOINC 18332-7]
<b>Methode, fabrikant:</b>	Fluorescence Enzyme Immuno Assay (FEIA) – ImmunoCAP Phadia 250, ThermoFisher Scientific

<b>Traceerbaarheid:</b>								
Gekalibreerd naar	Afkomstig uit de bijsluiter: Reference material: The ImmunoCAP Specific IgG Calibrators are traceable (via an unbroken chain of calibrations) to the International Reference Preparation (IRP) 67/86 for Human Serum Immunoglobulins A, G, and M from World Health Organisation (WHO). The calibrators are used for determination of Autoantigen specific IgG antibodies. The values are expressed in mgA/l, where A represents antigen-specific antibodies. By using a conversion factor included in the software, the values are expressed in IU/ml. The conversion factor has been determined by means of NIBSC research standard for anti-thyroid microsome serum code 66/387.							
<b>Referentie-interval of afkapgrenzen</b> (gebruik UCUM-eenheden)								
Herkomst referentiewaarden	Afkomstig van de firma: <table border="1" data-bbox="703 987 1348 1093"> <thead> <tr> <th>Assay</th> <th>Unit</th> <th>negative</th> </tr> </thead> <tbody> <tr> <td>ImmunoCAP Thyroid peroxidase IgG</td> <td>IU/mL</td> <td>&lt;60</td> </tr> </tbody> </table>		Assay	Unit	negative	ImmunoCAP Thyroid peroxidase IgG	IU/mL	<60
Assay	Unit	negative						
ImmunoCAP Thyroid peroxidase IgG	IU/mL	<60						
<b>Stabiliteit monster</b>								
Kamertemperatuur 4 °C -20 °C	Ja/Nee*, termijn: 2 dagen Ja/Nee*, termijn: 2 weken Ja/Nee*, termijn: 12 maanden							
<b>Detectielimieten</b> (gebruik UCUM-eenheden)								
LoD (Limit of Detection) LoQ (Limit of Quantitation)	N.v.t. bij autoantistofbepalingen							
<b>Imprecisie</b> (gebruik UCUM-eenheden)								
Concentratie (level invullen)	Level interne controle: 278 mg/L, CV (%): 8,8%							
<b>Meetbereik</b> (gebruik UCUM-eenheden)								
Meetbereik	33 – 3340 IU/mL							
<b>Extern QC programma</b>								
(Inter)nationaal extern QC programma Indien niet beschikbaar, alternatieve aanpak	SKML-HIM Schildklier antistoffen							

\* Doorhalen wat niet van toepassing is.

Ingevuld door: Caroline Roozendaal	Datum: 14-12-2018
------------------------------------	-------------------