

Performance data sheet

Testnaam Normetanefrine/kreatinine [molaire ratio] in urine

LOINC code 34348-3

Detection limit 0,1 umol/l

Imprecision

Within run:
Low 0,3 umol/l: CV% 1,25
Mid 3 umol/l: CV% 1,29
High 30 umol/l: CV% 0,93

Between run:
Low 0,3 umol/l: CV% 9,57
Mid 3 umol/l: CV% 6,11
High 30 umol/l: CV% 4,45

Measurement range 0,1 - 30 umol/l

Traceability

DL-normetanefrine:HCl MW: 219.7, Sigma N7127 / 100 mg
DL-normetanefrine:HCl-d3 MW: 222.7 Cambridge Isotopic Laboratories