

Clinical Laboratory Maxima Medisch Centrum

Assay performance data – Beta Trace Protein (BTC) in serum

Traceability:

- β -Trace Protein Internal protein reference preparation by Siemens (Behring)
- No international reference preparation is available

Detection limit:

- Serum (@ sample dilution 1:100) : 0,22 mg/L

Imprecision

- Serum in 1:100 sample dilution
At 0,22 mg/L : CV 12,9%
At 1.149 mg/L : CV 2,6%

Measurement range:

- Serum 0,22 (@1:100) - 35mg/L (@1:400)