

Performance data sheet

Testnaam	Melatonine [massa/volume] in speeksel
LOINC code	78992-5
Detection limit	<1,0 pg/ml
Imprecision	Low 6,836 pg/ml: CV% 7,7 Mid 13,274 pg/ml: CV% 8,4 High 20,029 pg/ml: CV% 6,5
Measurement range	<1,0 - >60,0 pg/ml
Traceability	Melatonine , MW 232,28 Sigma-Aldrich Cerilliant M-095, 1.0 mg/ml in 1 ml Methanol Melatonin-D7 solution MW 239,32 Sigma-Aldrich M-094 100 µg/mL in methanol, ampule of 1 mL